# waterpik®

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- ✓ Classic
- ✓ Implant Patient
- ✓ Periodontal Pocket
- ✓ Orthodontic Patient

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

(Carranza's Clinical Periodontology 12th Edition; Chapter 47; Page 510)

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
## Clinical Research Highlights

All Water Pik, Inc. clinical research is conducted by independent universities and research facilities and is published in peer reviewed journals. Recent research proves conclusively that the Waterpik® Water Flosser is more effective than string floss for removing plaque and reducing gingivitis and gingival bleeding.

### Waterpik® Water Flosser: Up to 29% More Effective for Plaque Removal than String Floss

Evaluation of the Plaque Removal Efficacy of a Water Flosser Compared to String Floss in Adults After a Single Use


Goyal CR, Lyle DM, Qaqish JG, Schuller R. et al. *J Clin Dent* 2013; 24(2):37-42. Study conducted at BioSci Research Canada, Ltd., Mississauga, Ontario, Canada.

### Waterpik® Water Flosser: Removes up to 99.9% of Plaque Biofilm After 3-Second Treatment

Biofilm Removal with a Dental Water Jet


Gorur A, Lyle DM, Schaudinn C, Costerton JW. *Compend Contin Ed Dent* 2009; 30 (Suppl 1):1-6. Study conducted at the University of Southern California School of Dentistry, USC Center for Biofilms, Los Angeles, California.

### Waterpik® Water Flosser: Up to 50% More Effective than String Floss for Reducing Gingivitis

Comparison of Irrigation to Floss as an Adjunct to Toothbrushing: Effect on Bleeding, Gingivitis and Supragingival Plaque


Barnes CM, Russell CM, Reinhardt RA et al. *J Clin Dent*, 2005; 16(3): 71-77. Study conducted at the University of Nebraska Medical Center, College of Dentistry, Lincoln, Nebraska.

## Waterpik® Water Flosser: Up to 80% More Effective than Air Floss for Reducing Gingivitis

Comparison of Two Power Interdental Cleaning Devices on the Reduction of Gingivitis


Sharma NC, Lyle DM, Qaqish JG, Schuller R. *J Clin Dent* 2012; 23(1): 22-26. Study conducted at BioSci Research Canada, Ltd., Mississauga, Ontario, Canada.

## Waterpik® Water Flosser: Up to 3X as Effective as String Floss for Orthodontic Patients

The Effect of a Dental Water Jet With Orthodontic Tip on Plaque and Bleeding in Adolescent Orthodontic Patients With Fixed Orthodontic Appliances


Sharma NC, Lyle DM, Qaqish JG et al. *J Ortho Dentofacial Orthop* 2008; 133(4): 565-571. Study conducted at BioSci Research Canada, Ltd., Mississauga, Ontario, Canada.

## Waterpik® Water Flosser: More Than 2X As Effective as String Floss for Implant Patients.

The Effect of a Water Flosser with Plaque Seeker Tip® on Gingival Bleeding for Implant Patients

Magnuson B, Harsono M, Stark PC, Lyle D, Kugel G, Perry R. Comparison of the effect of two interdental cleaning devices around implants on the reduction of bleeding: A 30-day randomized clinical trial. *Compend of Contin Ed in Dent* 2013; 34(Special Issue 8):2-7.

## Research Bibliography

1. Al-Mubarak S, Ciancio S, Aljada A, Awa H, Hamouida W, Ghanim H, Zambon J, Boardman T, Mohanty P, Ross C, Dandona P. *J Clin Periodontol* 2002; 29:295-300.
2. Aziz-Gandour IA, Newman HN. The effects of a simplified oral hygiene regime plus supragingival irrigation with chlorhexidine or metronidazole on chronic inflammatory periodontal disease. *J Clin Periodontol* 1986; 13:228-236.
3. Bakdash MB, et al. Daily compliance of chlorhexidine irrigation. Presented at IADR, Dublin, June 30, 1989.
4. Barnes CM, et al. Comparison of irrigation to floss as an adjunct to tooth brushing: Effect on bleeding, gingivitis and supragingival plaque. *J Clin Dent* 2005; 16(3):71-77.
5. Berger SA, et al. Bacteremia after the use of an oral irrigation device. *Annals of Int Med* 1974; 80:510-511.
6. Bhaskar S, et al. Water jet devices in dental practice. *J Periodontol* 1971; 42:658-664.
7. Bhaskar S, et al. Effect of high pressure water jet on oral mucosa of varied density. *J Periodontol* 1969; 40:593-598.
8. Boyd RL, et al. Comparison of a subgingivally placed cannula oral irrigator tip with a supragingivally placed standard irrigator tip. *J Clin Periodontol* 1992; 19:340-344.
9. Boyd RL, et al. Effect of self-administered daily irrigation with 0.02% SnF<sub>2</sub> on periodontal disease activity. *J Clin Periodontol* 1985; 12:420-431.
10. Brady JM, et al. Electronic microscopic study of the effect of water jet lavage on dental plaque. *J Dent Res* 1973; 52:1310-1313.
11. Braun RE, Ciancio SG. Subgingival delivery by an oral irrigating device. *J Periodontol* 1992; 63:469-472.
12. Brownstein CN, et al. Irrigation with chlorhexidine to resolve naturally occurring gingivitis. *J Clin Periodontol* 1990; 17:588-593.
13. Burch JG, et al. A two-month study of the effects of oral irrigation and automatic toothbrush use in an adult orthodontic population with fixed appliances. *Am J Orthod Dentofac Orthop* 1994; 106:121-126.
14. Cantor MT, Stahl SS. Interdental col tissue responses to the use of a water pressure cleansing device. *J Periodontol* 1969; 5:292-295.
15. Chaves ES, et al. Mechanism of irrigation effects on gingivitis. *J Periodontol* 1994; 65:1016-1021.
16. Ciancio SG, et al. Effect of a chemotherapeutic agent delivered by an oral irrigation device on plaque, gingivitis, and subgingival microflora. *J Periodontol* 1989; 60:310-315.
17. Cobb CM, et al. Ultrastructural examination of human periodontal pockets following the use of an oral irrigation device *in vivo*. *J Periodontol* 1982; 59:155-163.
18. Cutler CW, et al. Clinical benefits of oral irrigation for periodontitis are related to reduction of pro-inflammatory cytokine levels and plaque. *J Clin Periodontol* 2000; 27:134-143.
19. Dervanis JP, et al. Effects of a mouthwash in an irrigating device on accumulation and maturation of dental plaque. *J Periodontol* 1978; 49:81-84.
20. Drisko CL, et al. Comparison of dark-field microscopy and a flagella stain for monitoring the effect of a Water Pik on bacterial motility. *J Periodontol* 1987; 58:381-386.
21. Eakle WS, et al. Depth of penetration in to periodontal pockets with oral irrigation. *J Clin Periodontol* 1986; 13:39-44.
22. Felix JE, et al. Detection of bacteremia after the use of an oral irrigation device in subjects with periodontitis. *J Periodontol* 1971; 42:785-787.
23. Felo A, et al. Effects of subgingival chlorhexidine irrigation on peri-implant maintenance. *Am J Dent* 1997; 10:107-110.
24. Fine JB, et al. Short-term microbiological and clinical effects of subgingival irrigation with an antimicrobial mouthrinse. *J Periodontol* 1994; 65:30-36.
25. Flemmig TF, et al. Supragingival irrigation with 0.06% chlorhexidine in naturally occurring gingivitis. I. 6 month clinical observations. *J Periodontol* 1990; 61:112-117.
26. Flemmig TF, et al. Adjunctive supragingival irrigation with acetylsalicylic acid in periodontal supportive therapy. *J Clin Periodontol* 1995; 22:427-433.
27. Genovesi AM, et al. Periodontal maintenance following scaling and root planing. A randomized single-center study comparing minocycline treatment and daily oral irrigation with water. *Minerva Stomatol* 2013; 62(Suppl. 1 to NO. 12):1-9.
28. Gorur A, et al. Biofilm removal with a dental water jet. *Compend Contin Ed Dent* 2009; 30 (Suppl 1):1-6.
29. Goyal CR, et al. Evaluation of addition of a Water Flosser to Sonic Toothbrush. *J Clin Dent* 2012; 23(2):57-63.
30. Goyal CR, et al. Evaluation of the plaque removal efficacy of a water flosser compared to string floss in adults after a single use. *J Clin Dent* 2013; 24(2):37-42.
31. Hoover DR, Robinson HBG. The comparative effectiveness of a pulsating oral irrigator as an adjunct in maintaining oral health. *J Periodontol* 1971; 42:37-39.
32. Hugoson A. Effect of the Water Pik device on plaque accumulation and development of gingivitis. *J Clin Periodontol* 1978; 5:95-104.
33. Hurst JE, Madonia JV. The effect of an oral irrigating device on the oral hygiene of orthodontic patients. *JADA* 1970; 81:678-683.
34. Jolkovsky DL, et al. Clinical and microbiological effects of subgingival and gingival marginal irrigation with chlorhexidine gluconate. *J Periodontol* 1990; 61:663-669.
35. Krajewski J, et al. Evaluation of a water pressure cleaning device as an adjunct to periodontal treatment. *Periodontics* 1964; 2:76-78.
36. Kozam G. The effect of hydromassage on capillary strength. *NYSJDJ* 1973; 39:551-559.
37. Lainson PA, et al. A longitudinal study of pulsating water pressure cleansing devices. *J Periodontol* 1972; 43:444-446.
38. Larner JR, Greenstein G. Effect of calculus and irrigator tip design on depth of subgingival irrigation. *Int J Periodontics Res Dent* 1993; 13:289-297.
39. Lobene RR. The effect of a pulsed water pressure cleansing device on oral health. *J Periodontol* 1969; 40:667-670.
40. Magnuson B, et al. Comparison of the effect of two interdental cleaning devices around implants on the reduction of bleeding: A 30-day randomized clinical trial. *Compend Contin Ed Dent* 2013; 34(Special Issue 8):2-7.
41. Manhold JH, et al. Carbon penetration of gingival tissue by oral irrigating devices. *J Prev Dent* 1978; 5:3-6.
42. Newman MG, et al. Effectiveness of adjunctive irrigation in early periodontitis: Multi-center evaluation. *J Periodontol* 1994; 65:224-229.
43. Newman MG, et al. Irrigation with 0.06% chlorhexidine in naturally occurring gingivitis. II. 6-month microbiological observations. *J Periodontol* 1990; 61:427-433.
44. O'Leary TJ, et al. Possible penetration of crevicular tissue from oral hygiene procedures. I. *J Periodontol* 1970; 41:158-162.
45. Oshrain RL, et al. Oral irrigation devices: A clinical evaluation. *J Dent Hyg* 1987; 61:551-555.
46. Phelps-Sandall BA, Oxford SJ. Effectiveness of oral hygiene techniques on plaque and gingivitis in patients placed in intermaxillary fixation. *Oral Surg Oral Med Oral Pathol* 1983; 56:487-490.
47. Romans AR, App GR. Bacteremia, a result from oral irrigation in subjects with gingivitis. *J Periodontol* 1971; 42:757-760.
48. Rosema NAM, et al. The effect of different interdental cleaning devices on gingival bleeding. *J Int Acad Periodontol* 2011; 13(1):2-10.
49. Sander PC, et al. The effects of a simplified mechanical oral hygiene regime plus supragingival irrigation with chlorhexidine or metronidazole on subgingival plaque. *J Clin Periodontol* 1986; 13:237-242.
50. Schmidt, JC et al. Efficacy of various side-to-side toothbrushes for noncontact biofilm removal. *Clin Oral Invest* 2014; 18(3):793-800.
51. Schumacher HA. Experiences with the Water Pik irrigating device in multi-band treatment. *Quintessence Int* 1978; 9:47-49.
52. Selting WJ, et al. Water jet direction and periodontal pocket debridement. *J Periodontol* 1972; 43:569-572.
53. Sharma NC, et al. Comparison of two power interdental cleaning devices on plaque removal. *J Clin Dent* 2012; 23(1): 17-21.
54. Sharma NC, et al. Comparison of two power interdental cleaning devices on the reduction of gingivitis. *J Clin Dent* 2012; 23(1): 22-26.
55. Sharma NC, et al. Effect of a dental water jet with orthodontic tip on plaque and bleeding in adolescent patients with fixed orthodontic appliances. *Am J Ortho Dentofacial Orthop* 2008; 133(4):565-571.
56. Sharma NC, et al. Evaluation of the plaque removal efficacy of three power toothbrushes. *J Int Acad Periodontol* 2006; 8(3):83-88.
57. Tamimi GA, et al. Bacteremia study using a water irrigation device. *J Periodontol* 1969; 40:4-6.
58. Tempel TR, et al. Comparison of water irrigation and oral rinsing on clearance of soluble and particulate materials from the oral cavity. *J Periodontol* 1975; 46:391-396.
59. Waki MY, et al. Effects of subgingival irrigation on bacteremia following scaling and root planing. *J Periodontol* 1990; 61:405-411.
60. Walsh TF, et al. Clinical effects of pulsed oral irrigation with 0.2% chlorhexidine digluconate in patients with adult periodontitis. *J Clin Periodontol* 1992; 19:245-248.
61. Wolff LF, et al. The effect of professional and home subgingival irrigation with antimicrobial agents on gingivitis and early periodontitis. *J Dent Hyg* 1989; 63: 222-225, 241.
62. Wolff LF, et al. Phase contrast microscopic evaluation of subgingival plaque in combination with either conventional or antimicrobial home treatment of patients with periodontal inflammation. *J Periodontol Res* 1982; 17:537-540.

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